



MICROBIOLOGICAL HYGIENE CONTROL

BIOSOLUTE® microbiological culture media for hygiene monitoring of surfaces and air are available as ready-to-use plates, sedimentation plates and dehydrated culture media. Production is automated under clean room conditions in ISO-certified facilities. Strict quality controls in accordance with Ph. Eur. and ISO 11133 in the ISO 17025 certified test laboratory guarantee maximum product safety.

CONTACT PLATES, 55 MM, BLISTER, STORAGE AT RT

TOTAL BACTERIAL COUNT

Description	Standards	PK	Art. no.
Tryptic Soy Agar Agar (TSA)	Ph. Eur./USP, ISO	30 plates	9698.0030
Tryptic Soy Agar Agar TL	Ph. Eur./USP, ISO	30 plates	9691.0030
Tryptic Soy Agar Agar TLHTh	Ph. Eur./USP, ISO	30 plates	9321.0030
Plate Count Agar (PCA) TLHTh	DIN 10192, ISO 4833, ISO 8552, ISO 17410	30 plates	9553.0030
Plate Count Agar (PCA)	DIN 10192, ISO 4833, ISO 8552, ISO 17410	30 plates	9577.0030

FOR YEASTS AND MOLDS

Description	Standards	PK	Art. no.
DG 18 Agar	ISO 16000, ISO 21527	30 plates	9543.0030
DG 18 Metall-Agar	ISO 16000, ISO 21527	30 plates	9531.0030
Malt extract Agar	ISO 16000	30 plates	9982.0030
Malt extract Chloramphenicol Agar	ISO 16000	30 plates	8763.0030
Rose Bengal Chloramphenicol Agar (RBC Agar)		30 plates	9813.0030
Sabouraud Glucose Agar with TLHTh	Ph. Eur./USP, ISO 16212	30 plates	8161.0030
Sabouraud Dextrose Chloramphenicol Agar	Ph. Eur./USP, ISO 16212	30 plates	9225.0030

FOR COLIFORMS/ENTEROBACTERIA

Description	Standards	PK	Art. no.
MacConkey Agar No. 3	Ph. Eur./USP, ISO 21150	30 plates	9456.0030
VRB Agar (Crystal Violet Neutral Red Bile Lactose Agar)	ISO 4832, ISO 5541	30 plates	9666.0030
VRBD Agar (Violet Red Bile Dextrose Agar)	Ph. Eur./USP, ISO 21528	30 plates	9416.0030



Th. Geyer – we know lab

As a family business, we have been there for our customers for more than 130 years to support them in their demanding tasks in analytics, research and science.

Laboratories in industry, clinics, universities and colleges as well as private and public research institutions rely on our experience and ISO-certified quality when purchasing their laboratory supplies: from personal support to fast, secure delivery and important services, we reliably cover everything needed for complete laboratory equipment as a partner.

Th. Geyer Brands

Our Th. Geyer brand trio BIOSOLUTE®, CHEMSOLUTE® and LABSOLUTE® offers you attractive alternatives when it comes to life science products, laboratory chemicals and consumables.

BIOSOLUTE®

Products from microbiology, molecular biology, cell culture and histology specially tailored to your requirements in the life science laboratory.

Better orientation and easier selection for you in our constantly growing life science range. You can order all of these conveniently in our **webshop**.

FOR PSEUDOMONADS

Description	Standards	PK	Art. no.
Cetrimide Agar	Ph. Eur./USP, ISO 22717	30 plates	9246.0030

FOR STAPHYLOCOCCI

Description	Standards	PK	Art. no.
Baird Parker Agar	Ph. Eur./USP, ISO 5944, ISO 6888, ISO 22718	30 plates	9861.0030

PREPARED PLATES, 90 MM

TOTAL BACTERIAL COUNT

Description	Standards	PK	Art. no.
Tryptic Soy Agar (TSA)	Ph. Eur./USP, ISO	20 plates	9762.0020
Tryptic Soy Agar (TSA)	Ph. Eur./USP, ISO	120 plates	9762.0120

FOR YEASTS AND MOLDS

Description	Standards	PK	Art. no.
Rose Bengal Chloramphenicol Agar (RBC Agar)		20 plates	8666.0020
DG 18 Agar	ISO 16000, ISO 21527	20 plates	9526.0020
Malt extract Agar	ISO 16000	20 plates	9980.0020
Malt extract Chloramphenicol Agar	ISO 16000, ISO 16212	20 plates	9766.0020
Sabouraud 2% Dextrose Agar	Ph. Eur./USP, ISO 16212	20 plates	9784.0020
Sabouraud 2% Dextrose Chloramphenicol Gentamicin Cycloheximid Agar	Ph. Eur./USP, ISO 16212	20 plates	9174.0020
Sabouraud 4% Dextrose Chloramphenicol Agar	Ph. Eur./USP, ISO 16212	20 plates	9166.0020
Sabouraud Chloramphenicol-Gentamycin Agar	Ph. Eur./USP, ISO 16212	20 plates	8966.0020
Sabouraud Dextrose Agar	Ph. Eur./USP, ISO 16212	20 plates	9780.0020

FOR PSEUDOMONADS

Description	Standards	PK	Art. no.
Cetrimide Agar	Ph. Eur./USP, ISO 22717	20 plates	9248.0020

FOR COLIFORMS/ENTEROBACTERIA

Description	Standards	PK	Art. no.
MacConkey Agar Nr. 3	Ph. Eur./USP, ISO 21150	20 plates	9452.0020
VRB Agar (Crystal Violet Neutral Red Bile Lactose Agar)	ISO 4832, ISO 5541	20 plates	9875.0020
VRBD Agar (Violet Red Bile Dextrose Agar)	Ph. Eur./USP, ISO 21528	20 plates	9168.0020

FOR STAPHYLOCOCCI

Description	Standards	PK	Art. no.
Baird Parker Agar	Ph. Eur./USP, ISO 5944, ISO 6888, ISO 22718	20 plates	9746.0020
Baird Parker RPF Agar	Ph. Eur./USP, ISO 5944, ISO 6888, ISO 22718	20 plates	9748.0020

PREPARED PLATES, 90 MM, STORAGE AT ROOM TEMPERATURE

TOTAL BACTERIAL COUNT

Description	Standards	PK	Art. no.
Tryptic Soy Agar (TSA), double wrapping	Ph. Eur./USP, ISO	20 plates	8520.0020
R2A Agar RT, single wrapping	Ph. Eur./USP	20 plates	8527.0020

FOR YEASTS AND MOLDS

Description	Standards	PK	Art. no.
Sabouraud Dextrose Agar RT, double wrapping	Ph. Eur./USP, ISO 16212	20 plates	8521.0020

FOR COLIFORMS/ENTEROBACTERIA

Description	Standards	PK	Art. no.
Mac Conkey Agar No. 3 RT, single wrapping	Ph. Eur./USP, ISO 21150	20 plates	8524.0020
VRBD Agar RT, single wrapping	Ph. Eur./USP, ISO 21528	20 plates	8525.0020
XLD Agar Ph. Eur. RT, single wrapping	Ph. Eur./USP	20 plates	8526.0020

FOR PSEUDOMONADS

Description	Standards	PK	Art. no.
Cetrimide Agar RT, single wrapping	Ph. Eur./USP, ISO 22717	20 plates	8522.0020

FOR STAPHYLOCOCCI

Description	Standards	PK	Art. no.
Mannitol Salt Agar RT, single wrapping	Ph. Eur./USP, ISO 22718	20 plates	8523.0020

SETTLING PLATES, 90 MM

Description	Standards	PK	Art. no.
Tryptic Soy Agar (TSA) with TLHTh	Ph. Eur./USP, ISO	20 plates	9774.0020
Tryptic Soy Agar (TSA) TLHTh RT, triple wrapping, RD Box	Ph. Eur./USP, ISO	60 plates	9776.0060

DEHYDRATED CULTURE MEDIA AND SUPPLEMENTS

TOTAL BACTERIAL COUNT

Description	Standards	PK	Art. no.
Tryptic Soy Agar (TSA)	Ph. Eur./USP, ISO	500 g	9738.0500
Tryptic Soy Agar (TSA) with Polysorbate 80 and Lecithin		500 g	9775.0500
Plate Count Agar (PCA)	DIN 10192, ISO 4833, ISO 8552, ISO 17410	500 g	8425.0500

FOR YEASTS AND MOLDS

Description	Standards	PK	Art. no.
DG 18 Agar (Dichloran Glycerin Chloramphenicol Agar Base)	ISO 16000, ISO 21527	500 g	9685.0500
Malt extract agar	ISO 16000	500 g	8287.0500
Sabouraud 4 % Dextrose Agar	Ph. Eur./USP, ISO 16212	500 g	8135.0500
Chloramphenicol	Ph. Eur./USP	10 Vials	8778.0010

FOR COLIFORMS/ENTEROBACTERIA

Description	Standards	PK	Art. no.
MacConkey Agar No. 3	Ph. Eur./USP, ISO 21150	500 g	8796.0500
VRB Agar (Crystal Violet Neutral Red Bile Lactose Agar)	ISO 4832, ISO 5541	500 g	7883.0500
VRBD Agar (Violet Red Bile Dextrose Agar)	Ph. Eur./USP, ISO 21528	500 g	7836.0500

FOR PSEUDOMONADS

Description	Standards	PK	Art. no.
Cetrimide Agar (Pseudomonas Selective Agar)	Ph. Eur./USP, ISO 22717	500 g	9783.0500

FOR STAPHYLOCOCCI

Description	Standards	PK	Art. no.
Baird Parker Agar	Ph. Eur./USP, ISO 5944, ISO 6888, ISO 22718	500 g	9869.0500
Egg Yolk Tellurite Emulsion (20 %), sterile	ISO 5944, ISO 6888, ISO 22718	50 ml	9557.0050
Egg Yolk Tellurite Emulsion (20 %), sterile	ISO 5944, ISO 6888, ISO 22718	100 ml	9557.0100

DO YOU ALREADY KNOW ... ?

... the right consumables for hygiene control from LABSOLUTE®?

Our extensive range includes **wooden spatulas and cotton swabs** for sampling, **disposal bags and garbage bags** for safe disposal as well as **wiping and cleaning cloths** and **disposable gloves** for thorough cleaning.

CHEMSOLUTE® completes the range with reliable products for **surface disinfection** - ideally suited to your everyday laboratory hygiene measures.

You can find this and other articles in our webshop at www.thgeyer.com

For further information, please get in touch with your personal contact or send a short message to sales@thgeyer.de

