



SEMI-QUANTITATIVE TEST STRIPS

RELIABLE RESULTS IN JUST A FEW SECONDS

Our semi-quantitative test strips enable rapid yet reliable and precise analysis of water samples. The ready-to-use test systems require no calibration and are suitable for both laboratory use and on-site testing. Thanks to their ease of use and reproducible results, they support efficient routine analyses in industry, environmental monitoring, health-care, and research.

The practical snap-lid containers contain 100 test strips (*50 test strips).

Description	Measuring range	Art. no.
Ascorbic Acid 2000	0 – 50 – 100 – 200 – 300 – 500 – 700 – 1000 – 2000 mg/L	<u>2920.0001</u>
Chlorine 1000	0 – 50 – 100 – 250 – 500 – 1000 mg/L	<u>2971.0001</u>
Chlorine 20	0 – 0.5 – 1 – 3 – 5 – 10 – 20 mg/L	<u>2970.0001</u>
Total Hardness	3 – 4 – 7 – 14 – 21 °dH	<u>2915.0001</u>
3-in-1 Combination (pH, Chlorine, Alkalinity)*	pH 6.4 – 6.8 – 7.2 – 7.6 – 8.4, Chlorine 0 – 0.5 – 1 – 3 – 5 – 10 mg/L, Alkalinity 0 – 80 – 120 – 180 – 240 mg/L CaCO ₃	<u>3090.0001</u>
5-in-1 Combination (pH, Free and Total Chlorine, Alkalinity, Water Hardness)*	Total Hardness 0 – 100 – 250 – 500 – 1000 mg/L CaCO ₃ , Free Chlorine 0 – 0.5 – 1 – 3 – 5 – 10 mg/L Cl ₂ , Total Chlorine 0 – 1 – 3 – 5 – 10 mg/L Cl ₂ , Alkalinity 0 – 1 – 3 – 5 – 10 mg/L CaCO ₃ , pH 6.4 – 6.8 – 7.2 – 7.6 – 8.4	<u>3091.0001</u>
Nitrate 500	0 – 10 – 25 – 50 – 100 – 250 – 500 mg/L	<u>2935.0001</u>
Nitrite 25	0 – 0.5 – 1 – 5 – 10 – 25 mg/L	<u>2925.0001</u>
Nitrite 80	0 – 2 – 5 – 10 – 20 – 40 – 80 mg/L	<u>2930.0001</u>
Peracetic Acid 2000	0 – 500 – 1000 – 1500 – 2000 mg/L	<u>2947.0001</u>
Peracetic Acid 50	0 – 5 – 10 – 20 – 30 – 50 mg/L	<u>2945.0001</u>
Peracetic Acid 500	0 – 100 – 150 – 200 – 250 – 300 – 400 – 500 mg/L	<u>2946.0001</u>
Peroxide 100	0 – 1 – 3 – 10 – 30 – 100 mg/L	<u>2955.0001</u>
Peroxide 1000	0 – 100 – 200 – 400 – 600 – 800 – 1000 mg/L	<u>2941.0001</u>
Peroxide 25	0 – 0.5 – 2 – 5 – 10 – 25 mg/L	<u>2940.0001</u>
Phosphate 500	0 – 3 – 10 – 25 – 50 – 100 – 150 – 500 mg/L	<u>2950.0001</u>
QAC 1000	0 – 25 – 50 – 100 – 250 – 500 – 1000 mg/L	<u>2960.0001</u>

DO YOU ALREADY KNOW ... ?

... our LABSOLUTE® consumables and personal protective equipment?

You can find this and other articles in our webshop at www.thgeyer.com

For further information, please get in touch with your personal contact or send a short message to sales@thgeyer.de



Th. Geyer – We know Lab

As a family business, we have been there for our customers for more than 130 years to support them in their demanding tasks in analytics, research and science.

Laboratories in industry, clinics, universities and colleges as well as private and public research institutions rely on our experience and ISO-certified quality when purchasing their laboratory supplies: from personal support to fast, secure delivery and important services, we reliably cover everything needed for complete laboratory equipment as a partner.

Th. Geyer Brands

Our Th. Geyer brand trio BIOSOLUTE®, LABSOLUTE® and CHEMSOLUTE® offers you attractive alternatives when it comes to life science products, consumables and laboratory chemicals.

CHEMSOLUTE®

The perfect selection of chemicals for your demanding laboratory tasks: First-class laboratory chemicals, reagents and solutions – the perfect solution for your laboratory.

Laboratory chemicals from A-Z, measuring solutions, indicator solutions, chromatography reagents and much more. You can order all of these conveniently in our **webshop**.