



EFFICIENT ON-SITE MICROBIOLOGICAL MONITORING EASY SAMPLING, MEANINGFUL RESULTS

BIOSOLUTE® Dipslides are easy-to-use culture media for quickly and easily determining the microbiological contamination of liquids, surfaces, and process environments. By being directly inserted into liquids or swabbed on surfaces to be tested, they enable straightforward sampling without the need for laboratory work. Typical applications include the food industry, the beverage industry, water analysis, and pharmaceutical and cosmetics production. Additional benefits of BIOSOLUTE® Dipslides include a long shelf life of up to 7 months and convenient storage at room temperature (15 to 25°C).

ADVANTAGES

- **BIOSOLUTE® Dipslides:** Plastic paddles with different combinations of microbiological culture media
- **Easy to use:** Dip or press, incubate, evaluate
- **2 in 1:** Different combinations of selective agar on both sides
- **Fast results:** Reliable statements in a short time
- **High reproducibility:** Standardized media and clear growth patterns
- **Flexible use:** For water, cooling lubricants, process fluids and surfaces
- **Long shelf life:** Up to 7 months after date of manufacture
- **Storage at room temperature:** 15 °C – 25 °C

BIOSOLUTE® DIPSLIDES

Description	PK	Art. no.
Dipslides Total Count (PCA) / Disinfection Control (TSA with Tween(R) and Lecithin)	20 pieces	7890.0020
Dipslides Total Count (PCA) / Enterobacteria (VRBD)	20 pieces	7892.0020
Dipslides Total Count (PCA) / Coliforms (VRBL)	20 pieces	7894.0020
Dipslides Total Count (TSA with TTC) / Yeast and mold (Rose Bengal with Chloramphenicol)	20 pieces	7896.0020
Dipslides Total Count (Nutrient Agar) / Total Count (Nutrient Agar with TTC)	20 pieces	7898.0020



Th. Geyer – We know lab

As a family business, we have been there for our customers for more than 130 years to support them in their demanding tasks in analytics, research and science.

Laboratories in industry, clinics, universities and colleges as well as private and public research institutions rely on our experience and ISO-certified quality when purchasing their laboratory supplies: from personal support to fast, secure delivery and important services, we reliably cover everything needed for complete laboratory equipment as a partner.

Th. Geyer Brands

Our Th. Geyer brand trio BIOSOLUTE®, LABSOLUTE® and CHEMSOLUTE® offers you attractive alternatives when it comes to life science products, consumables and laboratory chemicals.

BIOSOLUTE®

Products from microbiology, molecular biology, cell culture and histology specially tailored to your requirements in the life science laboratory.

Better orientation and easier selection for you in our constantly growing life science range. You can order all of these conveniently in our **webshop**.