

LABSOLUTE® PCR STRIPS AND PLATES

COMPATIBILITY CHART



We confirm that LABSOLUTE® PCR strips and plates fit on following models*:

The information regarding the compatibility between the individual instrument and plate/strip is intended as a general guideline. Before making a purchase, you have the option to obtain samples of all plates/strips to ensure compatibility. We kindly ask you to verify the installed heat block and consult the instrument's manual for further details.

Skirt	PCR Strips									PCR Plates				
	-									non-skirted	semi-skirted			skirted
Number of tubes/wells	8	8	8	8	8	8	8	8	4	96	96	96	96	96
High/Low Profile (H/L)	H	H	H	L	L	H	L	H	-	H	H	H	L	L
Capacity (µl)	200	200	200	100	100	200	100	200	100	200	200	200	150	150
	7696510	7696511	7696514	7696512	7696515	7696516	7696517	7696518	7696519	7696535	7696540	7696541	7696545	7696530

ABI® / LIFE TECHNOLOGIES / THERMO FISHER SCIENTIFIC

Thermal Cyclers	96 well standard block	Veriti, Proflex, Simplicamp	x	x	x			x		x		x	x	x		
qPCR Cyclers	96 well standard block	GeneAmp® 2700/2720/9600/9700	x	x	x			x		x		x	x	x		
	96 well FAST block	7000, 7300, 7500, 7700, 7900 HT	x	x	x			x				x	x	x		
		QuantStudio™ 3/5/6/7/12K, ViiA7™	x	x	x			x				x	x	x		
		StepOne				x	x									
		StepOne Plus™				x	x		x							
		7500 FAST, 7900 HT FAST				x	x		x							
		QuantStudio™ 3/5/6/7/12K, ViiA7™				x	x		x							
Sequencer	96 well block	3100, 3130XL, 3500, 3500XL, 3730, 3730XL										x	x			

AGILENT / STRATAGENE

Thermal Cyclers	96 well block	Surecycler 8800	x	x	x							x				
qPCR Cyclers	96 well block	Robocycler Gradient				x	x		x			x				x
		AriaMx				x	x		x							x
		Mx3000P,™ Mx3005P™	x	x	x	x	x		x			x	x			
		Mx4000™	x	x	x							x	x			

ANALYTIK JENA / BIOMETRA

Thermal Cyclers	Strips only	TRIO, Tpersonal, T3 Thermocycler	x	x	x			x		x						
	96 well block	Flexcycler2, T1 Thermocycler, Tgradient, Tone, Tadvanced, TProfessional (Standard/Basic) Gradient/XL	x	x	x			x		x		x	x			
		Trobot 96, SpeedCycler2 (SP, SPR)	x	x	x	x	x	x	x			x	x			x
qPCR Cyclers	96 well block	qTOWER3/G/touch, Toptical	x	x	x			x				x	x			

LABSOLUTE® PCR STRIPS AND PLATES COMPATIBILITY CHART



We confirm that LABSOLUTE® PCR strips and plates fit on following models*:

The information regarding the compatibility between the individual instrument and plate/strip is intended as a general guideline. Before making a purchase, you have the option to obtain samples of all plates/strips to ensure compatibility. We kindly ask you to verify the installed heat block and consult the instrument's manual for further details.

<p>The information regarding the compatibility between the individual instrument and plate/strip is intended as a general guideline. Before making a purchase, you have the option to obtain samples of all plates/strips to ensure compatibility. We kindly ask you to verify the installed heat block and consult the instrument's manual for further details.</p>			Skirt	PCR Strips								PCR Plates					
				-								non-skirted	semi-skirted			skirted	
				8	8	8	8	8	8	8	8	4	96	96	96	96	96
				H	H	H	L	L	H	L	H	-	H	H	H	L	L
				200	200	200	100	100	200	100	200	100	200	200	200	150	150
				7696510	7696511	7696514	7696512	7696515	7696516	7696517	7696518	7696519	7696535	7696540	7696541	7696545	7696530
BIOER TECHNOLOGIES																	
Thermal Cyclers	Strips only	GeneQ	x	x	x			x		x							
	96 well block	Gene Touch 96	x	x	x			x		x		x	x				
BIO-RAD																	
Thermal Cyclers	Strips only	Genecycler	x	x	x			x		x							
	96 well block	C1000 Touch, S1000	x	x	x	x	x	x	x	x		x	x				x
		iCycler,™ MyCycler™, T100	x	x	x			x		x		x	x				
qPCR Cyclers	96 well block	CFX96 Touch, CFX96 Touch Deep Well, CFX connect				x	x		x								x
		MyiQ,™ iCycler™ IQ/IQ 4/IQ 5				x	x		x				x				
BIO-RAD MJ RESEARCH																	
Thermal Cyclers	Strips only	Mini Gradient	x	x	x			x		x							
	96 well block	PTC100™/200™/220™/221™/225™/240™	x	x	x	x	x	x	x	x		x	x				x
qPCR Cyclers	Strips only	MiniOpticon				x	x		x								
	96 well block	Chromo4™	x	x	x	x	x	x	x			x	x				x
		Opticon2™				x	x		x								x
CORBETT RESEARCH																	
Thermal Cyclers	96 well block	(Qiagen) Palm Cycler	x	x	x			x		x		x	x				
qPCR Cycler	Strips only	Rotor-Gene series									x						
CYTIVA / AMERSHAM																	
Sequencer	96 well block	MegaBACE™ 500, MegaBACE™ 1000 mark 2				x	x										x

LABSOLUTE® PCR STRIPS AND PLATES COMPATIBILITY CHART



We confirm that LABSOLUTE® PCR strips and plates fit on following models*:

The information regarding the compatibility between the individual instrument and plate/strip is intended as a general guideline. Before making a purchase, you have the option to obtain samples of all plates/strips to ensure compatibility. We kindly ask you to verify the installed heat block and consult the instrument's manual for further details.

			PCR Strips									PCR Plates			
Skirt			-									non-skirted	semi-skirted		skirted
Number of tubes/wells			8	8	8	8	8	8	8	8	4	96	96	96	96
High/Low Profile (H/L)			H	H	H	L	L	H	L	H	-	H	H	H	L
Capacity (µl)			200	200	200	100	100	200	100	200	100	200	200	200	150
			7696510	7696511	7696514	7696512	7696515	7696516	7696517	7696518	7696519	7696535	7696540	7696541	7696545
EPPENDORF															
Thermal Cyclers	96 well block	MasterCycler® ep/ep gradient/Pro/Pro S/nexus/nexus gradient/nexus SX1/nexus GSX1	x	x	x	x	x	x	x	x		x	x		x
		MasterCycler® nexus X2/GX2/GX2e/X2e	x	x	x	x	x	x	x	x		x			
qPCR Cycler	96 well block	Mastercycler™ ep realplex	x	x	x	x	x	x	x			x	x		x
PEQLAB / VWR															
Thermal Cyclers	Strips only	peqSTAR XS, peqSTAR 2X	x	x	x	x	x	x	x	x					
	96 well block	peqSTAR 96X	x	x	x	x	x	x	x	x		x	x		x
ROCHE															
qPCR Cyclers	96 well block	LC96, LC480				x	x								
	Strips only	Nano				x	x								
SENSOQUEST															
Thermal Cycler	96 well block	Labcycler	x	x	x	x	x	x	x	x		x	x	x	x
TAKARA															
Thermal Cycler	96 well block	Dice touch, Gradient	x	x	x			x		x		x	x		
TECHNE															
Thermal Cyclers	Strips only	3Prime, 3PrimeG, 3PrimeX	x	x	x			x		x					
	96 well block	Prime, PrimeG, Prime Elite, Prime Elite Satellite	x	x	x			x		x		x	x		
		PCRmax Alpha cycler 1/2/4	x	x	x			x		x		x	x		
		TC412, TC512, Genius, Genius Quad, Touchgene, Touchgene Gradient, Flexigene	x	x	x	x	x	x	x	x		x	x		x
qPCR Cycler	96 well block	Quantica				x	x		x						x