

# LABSOLUTE® PIPETTE TIPS COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

LABSOLUTE® Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
Single-channel pipette	0.1–2 µl 0.5–10 µl	1	x	x	x														
Single-channel pipette	2–20 µl	1				x		x	x	x									
Single-channel pipette	10–100 µl	1					x	x	x	x	x								
Single-channel pipette	20–200 µl	1						x	x	x	x								
Single-channel pipette	100–1000 µl	1											x	x					
Single-channel pipette	500–5000 µl	1													x				
Single-channel pipette	1–10 ml	1																x	
Multi-channel pipette	5–50 µl	8/12					x		x	x									
Multi-channel pipette	20–200 µl	8/12						x	x	x	x								
Multi-channel pipette	50–300 µl	8/12										x	x						

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

**LABSOLUTE® PIPETTE TIPS**  
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Sartorius/BioHit Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
mLine	0.1–3 µl 0.5–10 µl	1	x	x	x														
mLine	2–20 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
mLine	100–1000 µl	1											x	x					
mLine	0.5–10 µl	8/12	x	x	x														
mLine	5–100 µl	8/12				x	x	x	x	x	x								
mLine	30–300 µl	8/12										x	x						
Proline	0.1–2.5 µl 0.5–10 µl	1	x	x	x														
Proline	2–20 µl 5–50 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
Proline	100–1000 µl	1											x	x					
Proline	0.5–10 µl	8/12	x	x	x														
Proline	5–50 µl	8/12				x	x	x	x	x	x								
Proline	50–300 µl	8/12										x	x						

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

# LABSOLUTE® PIPETTE TIPS COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Brand Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
Transferpette	0.1–1 µl 0.5–10 µl	1																	
Transferpette	2–20 µl	1			x														
Transferpette	2–20 µl 5–50 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
Transferpette	100–1000 µl	1											x	x					
Transferpette	0.5–10 µl	8/12	x	x	x														
Transferpette	2–20 µl	8/12				x	x	x	x	x	x								
Transferpette	2.5–25 µl 5–50 µl 10–100 µl 20–200 µl	8/12				x	x	x		x	x								
Transferpette	30–300µl	8/12										x	x						
Transferpette	500–5000 µl	1															x		
Transferpette S	0.1–1 µl 0.1–2.5 µl 0.5–10 µl	1	x	x	x														
Transferpette S	2–20 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
Transferpette S	100–1000 µl	1											x	x					
Transferpette S	0.5–10 µl	8/12	x	x	x														
Transferpette S	5–50 µl 10–100 µl 20–200 µl	8/12				x	x	x	x	x	x								
Transferpette S	30–300 µl	8/12	x	x									x	x					
Transferpette S	500–5000 µl	1															x		

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

**LABSOLUTE® PIPETTE TIPS**  
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Eppendorf Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
Research Plus	0.1–2.5 µl 0.5–10 µl	1	x	x	x														
Research Plus	2–20 µl	1			x														
Research Plus	2–20 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
Research Plus	30–300 µl	1									x	x							
Research Plus	100–1000 µl	1											x	x					
Research Plus	0.5–10 µl	8/12	x	x	x														
Research Plus	10–100 µl	8/12				x	x	x	x	x	x								
Research Plus	30–300 µl	8/12										x	x						
Research Plus	500–5000	1														x			
Research Plus	1000–10000 µl	1																x	x
Reference	0.1–2.5 µl 0.5–10 µl	1	x	x	x														
Reference	2–20µl	1			x														
Reference	2–20 µl 10–100 µl 50–200 µl	1				x	x	x	x	x	x								
Reference	100–1000 µl	1											x	x					
Reference	500–5000 µl	1													x				
Reference	1000–10000 µl	1																x	x

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

**LABSOLUTE® PIPETTE TIPS**  
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Thermo/Finn Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
Finnpette 4500/ F1/F2	0.2 – 2 µl 0.5 – 5 µl	1	x		x			x	x	x		x							
Finnpette 4500/ F1/F2	0.5 – 5 µl 2 – 20 µl 5 – 50 µl 10 – 100 µl 20 – 200 µl	1				x	x	x	x	x	x	x							
Finnpette 4500/ F1/F2	100 – 1000 µl	1											x	x					
Finnpette 4500/ F1/F2	0.5 – 10 µl	8/12	x		x			x	x	x		x							
Finnpette 4500/ F1/F2	5 – 50 µl	8/12				x	x	x	x	x	x	x							
Finnpette 4500/ F1/F2	50 – 300 µl	8/12				x	x	x	x	x	x	x							
Finnpette 4500/ F1/F2	500 – 5000 µl	1															x		
Finnpette 4500/ F1/F2	1000 – 10000 µl	1																	x

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

**LABSOLUTE® PIPETTE TIPS**  
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Gilson Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
Pipetman Classic, G, L: P2	0.2–2 µl 1–10 µl	1	x	x	x														
Pipetman Classic, G, L: P20	2–20 µl 20–100 µl 50–200 µl	1				x	x	x	x	x	x								
Pipetman Classic, G, L: P1000	200–1000 µl	1											x	x					
Pipetman P10L	0.5–10 µl 2–20 µl 20–200 µl	8/12				x	x	x	x	x	x								
Pipetman P300L	20–300 µl	8/12										x	x						
Pipetman Classic, L, G	500–5000 µl	1														x			
Pipetman Classic, L, G	1000–10000 µl	1																	x

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

**LABSOLUTE® PIPETTE TIPS**  
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Nichiryo Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
Nichipet Premium, EX, EX Plus	0.1–2 µl	1	x	x	x														
Nichipet Premium, EX, EX Plus	2–20 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
Nichipet Premium, EX, EX Plus	100–1000 µl	1										x	x						
Nichipet 7000	0.5–10 µl	8/12	x	x	x														
Nichipet 7000	5–50 µl 40–200 µl	8/12				x	x	x	x	x	x								
Nichipet 7000	50–300 µl	8/12										x	x						

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

**LABSOLUTE® PIPETTE TIPS**  
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Rainin Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml		
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851
Pipet-Lite SL-2XLS	0.1–2 µl 0.5–10 µl	1	x	x	x														
Pipet-Lite SL-20XLS	2–20 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
Pipet-Lite SL-300XLS	20–300 µl	1										x	x						
Pipet-Lite SL-1000XLS Classic PR-2 Classic PR-20 Classic PR-1000	100–1000 µl	1											x	x					
Classic PR-2	0.1–2 µl 0.5–10 µl	1	x	x	x														
Classic PR-20	2–20 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x								
Classic PR-1000	100–1000 µl	1											x	x					
Classic PR-5000	500–5000 µl	1														x			
Classic PR-10000	5000–10000 µl	1																	x

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021



**LABSOLUTE® PIPETTE TIPS**  
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models\*:

Socorex Pipette	Volume range	Compatibility: Single-channel (1), multi-channel (8/12) pipette	10 µl		10 µl XL/ 20 µl	20 µl	100 µl	200 µl			300 µl		1,000 µl	1,250 µl	5 ml		10 ml			
			7695840, 7695860, 7695870, 7695890	7695880	7695841, 7695861, 7695871, 7695891, 7695881	7695882	7695883	7695842, 7695843	7695844, 7695845	7695862, 7695863, 7695872, 7695873, 7695892	7695884	7695846, 7695864, 7695874, 7695893	7695885	7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886	7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868	7695850	7695933	7695855, 7695856	7695851	7695934
Calibra® 822	0.2–2 µl 1–10 µl	1	x	x	x															
Calibra® 822	1–10 µl 2–20 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x									
Calibra® 822	100–1000 µl	1											x	x						
Calibra® 852	1–10 µl	8/12	x	x	x															
Calibra® 852	10–100 µl 20–200 µl	8/12				x	x	x	x	x	x									
Acura 825	0.1–2 µl 0.5–10 µl	1	x	x	x															
Acura 825	1–10 µl 2–20 µl 5–50 µl 10–100 µl 20–200 µl	1				x	x	x	x	x	x									
Acura 825	100–1000 µl	1											x	x						
Acura 855	0.5–10 µl	8/12	x	x	x															
Acura 855	5–50 µl 20–200 µl	8/12				x	x	x	x	x	x									
Acura 855	40–350 µl	8/12																		
Acura 835	500–5000 µl	1																	x	
Acura 835	1000-10000 µl	1																		x

\* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021