


LABSOLUTE® BECHERGLÄSER

low form, with spout

Properties:

- Made of borosilicate glass 3.3
- According to DIN 12331
- With functional spout and white, durable labelling field
- Dishwasher proof and autoclaveable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C

Technical drawing / picture:



The technical drawing is just an example. Technical drawings of all sizes are available on request.

Th. Geyer GmbH & Co. KG

Value table:

Item no.	V ml	Scale	D ml	Ø mm	h mm
7.690 000	5	-	-	22	30
7.690 001	10	-	-	26	35
7.690 002	25	10 – 20	10	34	50
7.690 003	50	10 – 25 – 40	15	42	60
7.690 004	100	25 – 50 – 75	25	50	70
7.690 005	150	50 – 75 – 100 – 125	25	60	80
7.690 006	250	50 – 100 – 150 – 200	50	70	95
7.690 007	400	100 – 200 – 300	50	80	110
7.690 008	600	100 – 200 – 300 – 400 – 500	100	90	125
7.690 009	800	200 – 300 – 400 – 500 – 600 – 700	100	100	135
7.690 010	1,000	200 – 400 – 600 – 800	100	105	145
7.690 011	2,000	500 – 750 – 1000 – 1250 – 1500	250	130	185
7.690 012	3,000	500 – 1000 – 1500 – 2000 – 2500	500	150	210
7.690 013	5,000	1000 – 2000 – 3000 – 4000	500	170	270

Description of the abbreviations in the value table:

Item no.	Item number
V	Nominal volume of the beaker in milliliter (ml)
Scale	Scaling that is printed on the beaker
D	Division of the scale in milliliter (ml)
Ø	Max. diameter of the beaker in millimeter (mm)
h	Max. height of the beaker in millimeter (mm)

Further information is part of the declaration of conformity.

Revision 1.0, Stand: 20.07.2017

Th. Geyer GmbH & Co. KG