


LABSOLUTE® BEAKERS

high form, with spout

Properties:

- Made of borosilicate glass 3.3
- According to DIN 12331
- With functional spout and white, durable labelling field
- Dishwasher proof and autoclaveable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C

Technical drawing / picture:



The technical drawing is just an example. Technical drawings of all sizes are available on request.

Th. Geyer GmbH & Co. KG

Value table:

Item no.	V ml	Scale	D ml	Ø mm	h mm
7.690 014	25	10 – 20	10	30	55
7.690 015	50	10 – 25 – 40	15	38	70
7.690 016	100	25 – 50 – 75	25	48	80
7.690 017	150	50 – 75 – 100 – 125	25	54	95
7.690 018	250	50 – 100 – 150 – 200	25	60	120
7.690 019	400	100 – 150 – 200 – 250 – 300	50	70	130
7.690 020	600	100 – 200 – 300 – 400 – 500	100	80	150
7.690 021	800	200 – 300 – 400 – 500 – 600 – 700	100	90	175
7.690 022	1,000	200 – 300 – 400 – 500 – 600 – 700 – 800	100	95	180
7.690 023	2,000	500 – 750 – 1000 – 1250 – 1500 – 1750	250	120	240
7.690 024	3,000	500 – 1000 – 1500 – 2000 – 2500	250	135	280

Description of the abbreviations in the value table:

Item no.	Item number
V	Nominal volume of the beaker in milliliter (ml)
Scale	Scaling that is printed on the beaker
D	Division of the scale in milliliter (ml)
Ø	Max. diameter of the beaker in millimeter (mm)
h	Max. height of the beaker in millimeter (mm)

Further information is part of the declaration of conformity.

Revision 1.0, Stand: 20.07.2017

Th. Geyer GmbH & Co. KG