

LABSOLUTE® AUTOMATIC BURETTES ACCORDING TO DR. SCHILLING

The automatic LABSOLUTE® burettes according to Dr. Schilling can be used for several applications in many areas like f. ex. electroplating, water treatment or sewage plants. The burettes of class B are justified to “Ex” and come with LDPE bottle and solid stand. The tolerances comply with DIN EN ISO 385. There is a Schellbach stripe on the backside of the burettes that makes it easier to read the exact filling level. The bottle and the burette are fixed together by the titrating armature.

Further features:

- precise dosing
- exact draining tip
- possibility of fine titration



Item no.	V _{Burette} ml	V _{Bottle} ml	Tolerance ± ml	Division ml	Pieces per pack
7 691 140	10	500	0.03	0.05	1
7 691 141	15	1000	0.04	0.1	1
7 691 142	25	1000	0.05	0.1	1
7 691 143	50	1000	0.08	0.1	1

Th. Geyer GmbH & Co. KG