

## **LABSOLUTE® CLOSURES FOR NARROW-MOUTH BOTTLES**

**fit to LABSOLUTE® items 7.696 161 to 7.696 165**

### **Properties:**

- Made of LDPE
- Food-safe according to Regulation (EC) No. 10/2011, 1282/2011 and 1183/2012, as well as Regulation (EC) No. 1935/2004
- Temperature resistance as usual for LDPE
- Chemical resistance as usual for LDPE

All closures are natural colored. Dropping caps have an additional red protection cap that preserves the tip and the liquid inside the bottle. The fine and long tip of the dropping caps allows to dose the liquid very precisely.

Standard screw caps close the bottles 100 % tight.

### **Value table standard screw caps:**

<b>Item no.</b>	<b>GL</b>
7.696 170	14
7.696 171	18
7.696 171	25

### **Value table dropping caps:**

<b>Item no.</b>	<b>GL</b>
7.696 174	14
7.696 175	18
7.696 176	25

### **Description of the abbreviations in the value table:**

Item no.      Item number  
GL            Thread size

**Revision 1.0, Stand: 20.03.2018**

**Th. Geyer GmbH & Co. KG**